



Heliogast Evolution

Heliogast Evolution is a Silicone adjustable gastric band. It is a one piece design with unlocking system and circular shape memory. It has an inflatable internal surface, totally atraumatic for the stomach. It has a mechanical closure system with hydraulic assistance.




Heliogast Advanced

Heliogast Advanced is a large Silicone adjustable gastric band. It is a one piece design with unlocking system and circular shape memory. It has an inflatable internal surface, totally atraumatic for the stomach. It has a mechanical closure system with hydraulic assistance.



Kit Heliogast

Kit Heliogast is a standard port replacement kit.



Kev1 Heliogast

Kev1 Heliogast is a smaller low profile port replacement kit which is used when patients have lost considerable weight and want to replace the regular port with a smaller one for aesthetic reasons. It has the same components as the regular port replacement kit.

