


Europclips

Titanium Atraumatic Haemostatic Clips, easy and convenient, soft loading system with sticky tape at the bottom of the cartridge. The chevron clip facilitates a safe "tip-first" closing action. Having a serrated inner surface and a heart shaped cross-section, the clip allows an atraumatic, but firm and permanent hold. Solid titanium gives a scatter-free CT image and allows for MIR compatibility.



Clip Appliers

Precise loading of the clips, firm jaw grip, safe closing, color coded, full range of sizes.



- AI-A-321530 Length 150mm straight
- AI-A-321532 Length 190mm straight
- AI-A-321534 Length 190mm right angled
- AI-A-321531 Length 190mm straight - super micro cardiac



- AI-A-321524 Length 150mm straight
- AI-A-321526 Length 190mm straight
- AI-A-321520 Length 190mm right angled
- AI-A-321522 Length 280mm straight
- AI-A-321528 Length 280mm right angled



- AI-A-321540 Length 190mm straight
- AI-A-321544 Length 190mm right angled
- AI-A-321542 Length 280mm straight
- AI-A-321546 Length 280mm right angled



- AI-A-321550 Length 190mm straight
- AI-A-321552 Length 190mm right angled
- AI-A-321554 Length 280mm straight
- AI-A-321556 Length 280mm right angled

Miscellaneous

Clip Appliers ... CONTINUED

Reusable Single Load Laparoscopic Clip Applicator, Clip Applicator for medium/large clips, rotatable.



Miscellaneous